



| Series | BT10 Series 10 µl Barrier Tip | 2040 Series 10 µl Ultra-Micro Tip | | |
|-------------------------------|--|--|---|---|
| Part Number | BT10.N | 2040.N | 2047.N 2047.NS | 2042.N 2042.NS |
| Graduation Marks | Graduation marks are indicated at both 0.5 µl, 2 µl and 10 µl locations – refer to product image for visuals | | | |
| Tip Composition | Neptune pipette tips are made of virgin polypropylene | | | |
| Tip Type | Natural Polypropylene | | | |
| Filter Material | High Density Polyethylene Filter | Non-filtered Products | | |
| Offered in Pre-Sterile Format | Yes | No | Yes | |
| Configuration | Racked | Bulk | ESP Reload* | Racked |
| Packaging Breakdown | 96 tips per rack 10 racks per pack 5 packs per case | 1000 tips per bag 20 bags per case | 96 tips per card 20 cards per pack 4 packs per case | 96 tips per rack 10 racks per pack 5 packs per case |
| Autoclavable | No | Autoclavable at 120 °C for 10-15 minutes at 15 PSI | | |
| Storage Conditions | Store in a clean, dry environment at room temperature 15-30 °C | | | |

BT10.N Series



2040.N Series



ESP Reload* Neptune's patented reload system
Environmentally Sustainable Pack (ESP) reduces plastic waste by up to 90%





Quality Control:

| | |
|----------------------------|---|
| Certificates of Compliance | Each lot undergoes stringent inspection and individual lot testing ensures Neptune products are certified RNase, DNase, DNA and Endotoxin-free. Visit www.neptunescientific.com to obtain a copy of a certificate of compliance for your Neptune product. |
| RNase/DNase | Products are washed in distilled water and concentrated via centrifugation. Samples are added to previously established nucleic acid standards, incubated for one hour at 37°C, and tested on a 2% gel using electrophoresis. Products must show no degradation of standards to pass. Test sensitivity is 10 ⁻⁷ Kunitz units/μl. |
| Nucleic Acid | Products are washed in distilled water and concentrated via centrifugation. Then, samples are added to protocol specified PCR reactions and thermal cycled for 50 cycles. A 2% agarose gel electrophoresis is used to examine experimental samples, positive controls, and negative controls. To pass, product samples must show no DNA amplification. Test sensitivity is 10 ng. |
| Endotoxin/Pyrogen | Products are tested for endotoxins by using the Limulus Amebocyte Lysate (LAL) gel assay according to FDA guidelines. Test sensitivity is 0.06 EU/ml. |
| Sterilization | Products are sterilized using electron beam irradiation. |
| Traceability | Each product contains a 5 digit lot number located on the rack, pack and case of each finished good. With Neptune's advanced manufacturing process all raw materials are able to be traced for maximum quality assurance. |

Advancements in Liquid Handling:

| | |
|---------------------|---|
| ESP Reload | Neptune's ESP (Environmentally Sustainable Pack) was the industry's first pipette reload system designed to minimize plastic waste by 90% and provide an environmentally friendly solution. |
| Aerosol Barrier Tip | Specifically engineered to reduce cross contamination. |

Pipettor Compatibility:

Biohit Proline™ 10 μl

Brand Transferpette S™ 10 μl

Capp™ 10 μl

CLP Beta-Pette™ 2 μl and 10 ul

CLP Poseidon™ 2 μl and 10 ul

CLP Poseidon Electronic™ 20 μl

Eppendorf Reference™ 2.5 μl

(works with 2040 series)

Eppendorf Reference™ 10 μl

Eppendorf Research™ 2.5 μl

(works with 2040 series)

Eppendorf Research™ 10 μl

Eppendorf Research Plus™ 2.5 μl

(works with 2040 series)

Eppendorf Research Plus™ 10 μl

Eppendorf Xplorer™ 10 μl

Finnpipette™ 10 μl and 50 μl

Finnpipette™ Electronic 10 μl

Finnpipette F1™ 10 μl

Gilson Pipetman™ P2 and P10

Hamilton™ 2 μl and 10 μl

Nichiryo Nichipet EX™ 10 μl

Nichiryo Oxford Benchmate™ 2 μl

Nichiryo Oxford Multimate™ 10 μl

Socorex Calibri 822™ 10 μl

VWR Ultra High Performance™

2 μl and 10 μl