

NEPTUNE®

Tools for Life Sciences

Optimus Rack Upgrade

Functional Tabs

- Simply squeeze and pull up for easy one handed snap card removal with reloads.

Locking Clasp

- Secures tips during transport.
- Clicks to ensure rack is securely closed.




Double Ball Hinge

- Upgraded hinge design withstands repeatable autoclaves.

Modified Hole Design

- Prevents rack lift up when removing tips.



Structural Support

- New SBS base footprint makes the rack more stable when attaching tips, eliminating the need for a Suretray.

Reusability

- The base is engineered with pyramidal wells and an SBS footprint making it reusable as a 96-well automation reservoir.

See what researchers and scientists have said about the Optimus Rack upgrades

"It is great! Very user-friendly."

"The hinged design works well with multichannel pipets and the design of the box keeps the tips from all falling out when you place them on your pipette."

"I love that it is eco-friendly! Helps lab be greener, BIG positive!!"

"The locking clasp is very sturdy, superior to other brands."

www.neptunescientific.com/optimus