

# PCR PLATE COMPATIBILITY CHART

Plate PN		3730.X	3438.X	3742.X	3732.X
Standard well (Overall height approx. 21 mm, max capacity > 350 µl)		●		●	
Low well (Overall height approx. 16 mm, max capacity > 200 µl)			●		●
<b>Applied Biosystems</b>					
Thermal Cyclers	2700	●		●	
	9600	●		●	
	9700	●		●	
	9800 "Fast"				
	Veriti 0.1mL				
	Veriti 0.1mL				
	Veriti 384				
"Real Time" Thermal Cyclers	5700	●		●	
	PRISM 7000	●		●	
	7300	●		●	
	7500	●		●	
	7500 "Fast"				
	7700	●		●	
	7900HT	●		●	
	Step One				
	Step One plus				
Sequencers	PRISM 310	●		●	
	PRISM 3100	●		●	
	3130 (XL)	●		●	
	3700	●		●	
	PRISM 3730 (XL)	●		●	
<b>Amersham</b>					
Sequencers	MegaBACE 500				●
	MegaBACE 1000				●
<b>Beckman</b>					
Sequencers	CEQ	●			
<b>Biometra</b>					
Thermal Cyclers	Uno	●	●	●	●
	Uno II	●	●	●	●
	T1 Thermal Cycler	●	●	●	●
	Tgradient	●	●	●	●
	Trobot	●	●	●	●
<b>Bio-Rad/MJ Research</b>					
Thermal Cyclers	Gene Cycler	●		●	
	PTC-100	●	●		●
	PTC-200	●	●		●
	PTC-225 Tetrad	●	●		●
	Dyad/Dyad Disciple	●	●		●
	iCycler	●			●
	Mycycler	●			
"Real Time" Thermal Cyclers	Opticon/Opticon2	●	●		●
	Chromo-4		●		●
	iCycler	●			●
	MyiQ	●			●
	iQ5	●			●
	CFX96				●
	CFX384				
Sequencers	BaseStation				●
<b>Corbett Research (Qiagen)</b>					
Thermal Cyclers	Palm Cycler 96		●		●
	Palm Cycler 384				
	RotorGene				

PCR Plate Compatibility Chart



Plate Part Number		3730.X	3438.X	3742.X	3732.X
<b>Eppendorf</b>					
Thermal Cyclers	Mastercycler	●	●		●
	Mastercycler Gradient	●		●	
	Mastercycler ep	●	●		●
	Mastercycler M384				
"Real Time" Thermal Cyclers	Mastercycler ep Realplex				●
<b>Ericom</b>					
"Real Time" Thermal Cyclers	Power Block I	●	●	●	
	Deltacycler I	●	●		
	Deltacycler II	●	●	●	
	Single Block	●	●		
	Twin Block	●	●		
<b>Esco</b>					
Thermal Cyclers	Swift	●			
<b>G-Storm</b>					
Thermal Cyclers	GS1	●	●	●	
	GS4	●	●	●	
	GSX	●	●	●	
	GSXs	●	●	●	
<b>MWG</b>					
Thermal Cyclers	Primus 96	●	●	●	●
	Primus 384				
<b>Stratagene</b>					
Thermal Cyclers	Robocycler 96	●			
	Robocycler Gradient	●	●	●	●
"Real Time" Thermal Cyclers	Mx4000	●	●	●	
	Mx3000	●	●	●	
<b>TaKaRa</b>					
Thermal Cyclers	TP240				●
	TP3000	●	●	●	●
<b>Techne</b>					
Thermal Cyclers	Touchgene	●	●	●	
	Cyclogene	●	●	●	
	Genius	●	●	●	
	Genius Quad	●	●	●	
	Genius (TC412)	●	●	●	●
	Flexigene	●	●	●	●
	Touchgene X	●	●		●
	Touchgene Gradient (TC512)	●	●	●	●
"Real Time" Thermal Cyclers	Quantica		●		●
<b>Thermo Hybrid</b>					
Thermal Cyclers	PCR Sprint	●	●	●	●
	MBS Satellite System	●	●	●	●
	Px2	●	●	●	●
	PxE	●	●	●	●
	PCR Express	●	●	●	●
	Omni-E	●	●	●	●
	Touchdown	●	●	●	●
OmniGene	●	●	●	●	
<b>Transgenomic</b>					
Sequencers	WAVE System				●