

REAGENT RESERVOIRS

White Box Disposable Reagent Reservoirs

White Box Disposable Reagent Reservoirs offer precision molded design, with the convenience and economy of a disposable. These disposable reservoirs are designed for use with multi-channel pipetting applications and can accommodate up to 12-channel pipettes. Angled sides and a V-shaped bottom allow maximum recovery of reagents. All reservoirs feature internal graduated volume lines and corner pour-spouts. Reservoirs are available in several packaging configurations for laboratory convenience.


As with all products in the White Box line, disposable reagent reservoirs are produced in an ISO 9001-2008 certified facility and each lot of product undergoes rigorous QC inspection, including laboratory analysis to certify against the presence of RNase, DNase, Endotoxin (pyrogens), Nucleic Acids and PCR inhibitors.



25ml Reagent Reservoir



**25ml Reagent Reservoir
3-Divided Wells**



50ml Reagent Reservoir



100ml Reagent Reservoir

WHITE BOX DISPOSABLE REAGENT RESERVOIRS

Reagent Reservoirs in Packs of Five

| Part Number | Description | Packaging |
|---------------|---|-----------------------------------|
| RR-25-5-PS | 25ml Reagent Reservoir, Polystyrene, Pre-Sterile | 5 Reservoirs/Bag; 20 Bags/Case |
| RR-25-5-PS-D3 | 25ml Reagent Reservoir, 3-Divided Wells, Polystyrene, Pre-Sterile | 5 Reservoirs/Bag; 20 Bags/Case |
| RR-50-5-PS | 50ml Reagent Reservoir, Polystyrene, Pre-Sterile | 5 Reservoirs/Bag; 20 Bags/Case |
| RR-100-5-PS | 100ml Reagent Reservoir, Polystyrene, Pre-Sterile | 5 Reservoirs/Bag; 20 Bags/Case |

Bulk Packaged Reagent Reservoirs

| Part Number | Description | Packaging |
|---------------|---|---------------------|
| RR-25-B-NS | 25ml Reagent Reservoir, Polystyrene, Non-Sterile | 100 Reservoirs/Case |
| RR-25-B-NS-D3 | 25ml Reagent Reservoir, 3-Divided Wells, Polystyrene, Non-Sterile | 100 Reservoirs/Case |
| RR-50-B-NS | 50ml Reagent Reservoir, Polystyrene, Non-Sterile | 100 Reservoirs/Case |
| RR-100-B-NS | 100ml Reagent Reservoir, Polystyrene, Non-Sterile | 100 Reservoirs/Case |

Individually Wrapped Reagent Reservoirs

| Part Number | Description | Packaging |
|----------------|---|----------------------------------|
| RR-25-IW-PS | 25ml Reagent Reservoir, Polystyrene, Pre-Sterile | 1 Reservoir/Bag; 50 Bags/Case |
| RR-25-IW-PS-D3 | 25ml Reagent Reservoir, 3-Divided Wells, Polystyrene, Pre-Sterile | 1 Reservoir/Bag; 50 Bags/Case |
| RR-50-IW-PS | 50ml Reagent Reservoir, Polystyrene, Pre-Sterile | 1 Reservoir/Bag; 50 Bags/Case |
| RR-100-IW-PS | 100ml Reagent Reservoir, Polystyrene, Pre-Sterile | 1 Reservoir/Bag; 50 Bags/Case |

A BUILDING BLOCK FOR YOUR BRANDED PORTFOLIO...

Interested in offering liquid handling products under your own corporate brand, but lacking the volume necessary for private label? Our White Box portfolio is the perfect solution.

Ideal for supplementing assay kits or expanding your product line into new market segments, White Box products are a cost effective solution and include the most widely used laboratory consumables required in life science and clinical research.

Manufactured under stringent controls and of the highest quality, products in the White Box line are packaged and shipped in discrete white boxes, ready for you to customize with your own branded labels.



Commercial Benefits of the White Box Product Line include:

- Sold in case quantities
- No purchasing minimums
- Maintained as standard inventoried items
- Reduced lead time
- High quality unit boxes and packaging

WHITE BOX AVAILABLE PRODUCTS



**Universal
Pipette Tips**



**Pipette Tip
Reload System**



**Disposable
Reagent Reservoirs**

Ordering and Technical Information:

White Box products are supported by a dedicated team of professionals with over 25 years of industry experience. Please call or email Customer Service to place an order or submit technical questions.

Phone: 858.875.7696 or 888.330.2396 (6am to 5pm USA PST)
Fax: 858.875.5496
Email: csr@biotix.com

Manufactured by:

