

## BULK PACKAGING

---

# White Box Universal Pipette Tips

Among the many options in the White Box product line, universal pipette tips are available in a bulk-packaged format. As an economical and flexible offering, pipette tips are packaged in re-sealable plastic bags to help avoid contamination over the course of use. To support broad research applications, bulk packaged pipette tips are available in both non-filtered and filtered versions. The filtered tips utilize premium, high-density polyethylene filters specifically designed to prevent cross contamination in PCR and other high sensitivity assays.









As with all products in the White Box line, bulk packaged universal pipette tips are produced in an ISO 9001-2008 certified facility and each lot of product undergoes rigorous QC inspection, including laboratory analysis to certify against the presence of RNase, DNase, Endotoxin (pyrogens), Nucleic Acids and PCR inhibitors. White Box universal pipette tips are designed to be compatible with most of the industry's leading pipettor brands. To download the pipette tip compatibility chart, please visit [neptunescientific.com/services/whitebox](http://neptunescientific.com/services/whitebox).

Phone: 858.875.7696 or 888.330.2396 (6am to 5pm USA PST)  
Fax: 858.875.5496  
Email: [csr@biotix.com](mailto:csr@biotix.com)

Manufactured by:   
Fluid Innovation







# WHITE BOX UNIVERSAL PIPETTE TIPS: BULK PACKAGING

## Universal Pipette Tips, Non-Filtered

Part Number	Description	Packaging
<b>P10B</b>  <small>Actual Size</small>	10µl Micro Tip, Clear	1,000 Tips/Pack, 20 Packs/Case
<b>P200B</b>  <small>Actual Size</small>	200µl Pipette Tip, Clear	1,000 Tips/Pack, 10 Packs/Case
<b>P200BL</b>  <small>Actual Size</small>	200µl Low Retention Pipette Tip, Clear	1,000 Tips/Pack, 10 Packs/Case
<b>P2700-200</b>  <small>Actual Size</small>	200µl Pipette Tip, Yellow	1,000 Tips/Pack, 10 Packs/Case
<b>P1000B</b>  <small>Tip image is 85% of actual size</small>	1000µl Pipette Tip, Clear	1,000 Tips/Pack, 4 Packs/Case
<b>P2850-1MBB</b>  <small>Tip image has been reduced to 84% of actual size</small>	1000µl Pipette Tip, Blue	1,000 Tips/Pack, 4 Packs/Case

# WHITE BOX UNIVERSAL PIPETTE TIPS: BULK PACKAGING

## Universal Pipette Tips, Filtered

Part Number	Description	Packaging
<b>F10B</b>  <small>Actual Size</small>	10µl Filtered Micro Tip, Clear	1,000 Tips/Pack, 10 Packs/Case
<b>F10XB</b>  <small>Actual Size</small>	10µl Filtered Extended Length Tip, Clear	1,000 Tips/Pack, 10 Packs/Case
<b>F20B</b>  <small>Actual Size</small>	20µl Filtered Pipette Tip, Clear	1,000 Tips/Pack, 10 Packs/Case
<b>FT100-B.X</b>  <small>Actual Size</small>	100µl Filtered Pipette Tip, Clear	1,000 Tips/Pack, 10 Packs/Case
<b>F200B</b>  <small>Actual Size</small>	200µl Filtered Pipette Tip, Clear	1,000 Tips/Pack, 10 Packs/Case
<b>FT1000-B.X</b>  <small>Tip image has been reduced to 85% of actual size</small>	1000µl Filtered Pipette Tip, Clear	1,000 Tips/Pack, 4 Packs/Case

# A BUILDING BLOCK FOR YOUR BRANDED PORTFOLIO...

**Interested in offering liquid handling products under your own corporate brand, but lacking the volume necessary for private label? Our White Box portfolio is the perfect solution.**

Ideal for supplementing assay kits or expanding your product line into new market segments, White Box products are a cost effective solution and include the most widely used laboratory consumables required in life science and clinical research.

Manufactured under stringent controls and of the highest quality, products in the White Box line are packaged and shipped in discrete white boxes, ready for you to customize with your own branded labels.



## **Commercial Benefits of the White Box Product Line include:**

- Sold in case quantities
- No purchasing minimums
- Maintained as standard inventoried items
- Reduced lead time
- High quality unit boxes and packaging

## **WHITE BOX AVAILABLE PRODUCTS**



**Universal  
Pipette Tips**



**Pipette Tip  
Reload System**



**Disposable  
Reagent Reservoirs**

## **Ordering and Technical Information:**

White Box products are supported by a dedicated team of professionals with over 25 years of industry experience. Please call or email Customer Service to place an order or submit technical questions.

Phone: 858.875.7696 or 888.330.2396 (6am to 5pm USA PST)  
Fax: 858.875.5496  
Email: [csr@biotix.com](mailto:csr@biotix.com)

Manufactured by:

