

TIP RELOAD SYSTEM

White Box Universal Pipette Tip Reload System



Economical and environmentally friendly, the White Box Reload System reduces plastic and packaging waste by up to 90%. Tips in the reload format provide a low cost alternative compared to individually racked product, while saving time not having to manually load bulk tips. This transfer system allows you to reload your empty trays with new tips in a single movement.

Key features of the White Box Tip Reload System include:

- Requires 57% less storage space
- Reuse of existing tip trays
- Easy, single-handed reload without transfer tools

As with all products in the White Box line, the tip reload systems are produced in an ISO 9001-2008 certified facility and each lot of product undergoes rigorous QC inspection, including laboratory analysis to certify against the presence of RNase, DNase, Endotoxin (pyrogens), Nucleic Acids and PCR inhibitors. White Box universal pipette tips are designed to be compatible with most of the industry's leading pipettor brands. To download the pipette tip compatibility chart, please visit neptunescientific.com/services/whitebox.






Phone: 858.875.7696 or 888.330.2396 (6am to 5pm USA PST)
Fax: 858.875.5496
Email: csr@biotix.com

Manufactured by:

Biotix
Fluid Innovation

WHITE BOX UNIVERSAL PIPETTE TIP RELOAD SYSTEM

Universal Pipette Tips Available in Reload System

Part Number	Description	Packaging
10UL-RELOAD  <small>Actual Size</small>	10µl Micro Tip, Clear, Low Retention	20 Cards of 96 Tips/Pack, 4 Packs/Case
10XL-RELOAD  <small>Actual Size</small>	10µl Extended Length Tip, Clear, Low Retention	10 Cards of 96 Tips/Pack, 10 Packs/Case
200UL-RELOAD  <small>Actual Size</small>	200µl Pipette Tip, Clear, Low Retention	10 Cards of 96 Tips/Pack, 10 Packs/Case
1000UL-RELOAD  <small>Tip image is 85% of actual size</small>	1000µl Pipette Tip, Clear, Low Retention	10 Cards of 96 Tips/Pack, 5 Packs/Case
1250UL-RELOAD  <small>Tip image is 72% of actual size</small>	1250µl Extended Length Pipette Tip, Clear, Low Retention	10 Cards of 96 Tips/Pack, 5 Packs/Case

Empty Tip Trays for White Box Reload System

Part Number	Description	Packaging
EMTRSM LN4	Empty Trays for 10µl and 200µl Tip Reloads	50 Trays per Case
EMTRLARN4	Empty Trays for 1000µl and 1250µl Tip Reloads	40 Trays per Case

A BUILDING BLOCK FOR YOUR BRANDED PORTFOLIO...

Interested in offering liquid handling products under your own corporate brand, but lacking the volume necessary for private label? Our White Box portfolio is the perfect solution.

Ideal for supplementing assay kits or expanding your product line into new market segments, White Box products are a cost effective solution and include the most widely used laboratory consumables required in life science and clinical research.

Manufactured under stringent controls and of the highest quality, products in the White Box line are packaged and shipped in discrete white boxes, ready for you to customize with your own branded labels.



Commercial Benefits of the White Box Product Line include:

- Sold in case quantities
- No purchasing minimums
- Maintained as standard inventoried items
- Reduced lead time
- High quality unit boxes and packaging

WHITE BOX AVAILABLE PRODUCTS



**Universal
Pipette Tips**



**Pipette Tip
Reload System**



**Disposable
Reagent Reservoirs**

Ordering and Technical Information:

White Box products are supported by a dedicated team of professionals with over 25 years of industry experience. Please call or email Customer Service to place an order or submit technical questions.

Phone: 858.875.7696 or 888.330.2396 (6am to 5pm USA PST)
Fax: 858.875.5496
Email: csr@biotix.com

Manufactured by:

